

QuikCoag PT-HS

The **QuikCoag PT-HS** reagent is an in vitro diagnostic reagent intended for use in performing the one-stage prothrombin time (PT) test and assays which are based on a modified prothrombin time.

- ◆ Rabbit brain based reagent
- ◆ Convenient lyophilized format
- ◆ Reagent can be used on any open system
- ◆ Available in 2.0 mL, 4.0 mL and 10.0 mL fill volumes
- ◆ Reconstituted reagent is stable for five days when stored in the original container at 2 to 8°C
- ◆ Three year expiry from the date of manufacture
- ◆ Target ISI value range is 0.90 – 1.50, with a specific ISI assigned to each lot
- ◆ CE Marked, FDA 510(k) cleared and Health Canada Registered
- ◆ Competitively priced
- ◆ Available worldwide through a network of carefully selected distributors
- ◆ Unique Device Identifiers with scannable and human readable bar codes on box and vial labels
- ◆ OEM and customized options available

Related Products

- ◆ QuikCoag Control Level 1 Normal
- ◆ QuikCoag Control Level 2 Low Abnormal
- ◆ QuikCoag Control Level 3 High Abnormal
- ◆ QuikCoag MP Control Level 1 Normal
- ◆ QuikCoag MP Control Level 2 Abnormal



biomedicadiagnostics.com

94 Wentworth Road, PO Box 1030, Windsor, NS Canada B0N 2T0 T 902 798 5105 F 902 798 1025