

## Medicinal Freezer

## MED110MF

The MED110MF chest freezer is specifically designed for modern hospital and pharmacy facilities to preserve medical products to tight temperature tolerances. A digital controller is employed to maintain freezer set temperatures between 0°C and -40°C whilst having an audible alarm function should set temperature parameters be exceeded. Whilst the exterior casing is made from white Chromadek steel, the interior finish is in stainless steel for an enhanced hygiene function. The MED110MF has the added advantage to conform to the latest environmental 'green' standards for gases and electrical consumption. This product is backed by the reputable and reliable Minus 40 service warranty and customers can select maintenance contracts to ensure extended longevity of their Minus 40 Fridge.



### SPECIFICATIONS

Volume	110L
Temperature Range	0°C to - 40°C
Suitable for maximum ambient temperature	35°C
Voltage	AC 220V/50Hz
Refrigerant	R507a

### FEATURES

- The outer casing is produced with white pre-painted 0,7mm Chromadek.
- The inner casing is produced with stainless steel.
- This unit's insulation is foamed in place by Polyurethane 40kg/m<sup>2</sup> density with a minimum of 80mm thickness.
- Control of the refrigerator temperature is by means of a digital thermostat.
- Alarm indication can be installed for temperatures above 0°C.
- An environmentally friendly refrigerator that is CFC free.
- This unit consumes minimum power.
- Chest lid is lockable door for product security.
- An operating and installation manual is provided.

### EXTERNAL DIMENSIONS

Height mm	900
Width mm	915
Depth mm	560
Weight kg	80

### INTERNAL DIMENSIONS

Height mm	400
Width mm	750
Depth mm	360

