



Vaccine Freezer Comparison

	MED60MF	MED110MF	MED230MF	MED540MF	MB400LF	BB470LF
VOLUME	60L	110L	230L	540L	540L	470L
TEMPERATURE RANGE	0°C to - 40°C	0°C to - 40°C	-10°C to - 40°C	0°C to - 25°C	0°C to - 25°C	0°C to - 40°C
AMBIENT TEMPERATURE	35°C	35°C	35°C	35°C	35°C	35°C
VOLTAGE	AC220V/50Hz	AC220V/50Hz	AC220V/50Hz	AC220V/50Hz	AC220V/50Hz	AC220V/50Hz
REFRIGERANT	R507a	R507a	R507a	R404a	R404a	R507a