



Minus 40 Vaccine Fridge / Freezers

| | B130/10/12 | B130/10/220H | B200/40/22H |
|--|-------------------|---------------------|--------------------|
| FRIDGE COMPARTMENT VOLUME | 130L | 130L | 200L |
| FREEZER COMPARTMENT VOLUME | 30L | 30L | 40L |
| FRIDGE TEMPERATURE RANGE | 2°C to 8°C | 2°C to 8°C | 2°C to 8°C |
| FREEZER TEMPERATURE RANGE | -.12°C to -10°C | -.15°C to -10°C | -.15°C to -10°C |
| HOLD-OVER TIME OF + 10°C (32°C AMB) | > 6HRS | > 48HRS | > 55HRS |
| SUITABLE FOR MAXIMUM AMBIENT TEMPERATURE | 43°C | 43°C | 43°C |
| ICE PACK FREEZING PER KG / 24HRS | 4.8KG | 4.8KG | 8KG |
| VOLTAGE / FREQUENCY | 12V or 24V | 220V AV . 50HZ | 220V AV . 50HZ |
| ELECTRICITY CONSUMPTION WITHOUT ICE PACK FREEZING | 0.66kW / 24Hrs | | |
| REFRIGERANT | R134a | R134a | R134a |
| CLIMATE CLASS | T | T | T |
| DEFROST CYCLE | Manual | Manual | Manual |