

SCIENTIFIC MANUFACTURING DIGITAL OVENS**Model: 275**

Voltage V/hz	230/50
Power W	1000
Shipping Size WxHxD (mm)	810 x 550 x 670
Mass (kg)	38
Internal dimensions WxHxD (mm)	385 x 375 x 370
External dimensions WxHxD (mm)	651 x 510 x 600
Internal Volume (l)	40
Temp. variation across chamber ($\pm^{\circ}\text{C}$)	1.5
Temp. variation at a point ($\pm^{\circ}\text{C}$)	0.3
Number of Shelves	2
Number of Shelf Positions	5
Maximum load per shelf (kg)	26

Model: 277

Voltage V/hz	230/50
Power W	1500
Shipping Size WxHxD (mm)	970 x 950 x 840
Mass (kg)	75
Internal dimensions WxHxD (mm)	535 x 580 x 535
External dimensions WxHxD (mm)	801 x 720 x 770
Internal Volume (l)	160
Temp. variation across chamber ($\pm^{\circ}\text{C}$)	2
Temp. variation at a point ($\pm^{\circ}\text{C}$)	0.3
Number of Shelves	2
Number of Shelf Positions	8
Maximum load per shelf (kg)	26

Model: 278D

Voltage V/hz	230/50
Power W	2000
Shipping Size WxHxD (mm)	820 x 1090 x 920
Mass (kg)	88
Internal dimensions WxHxD (mm)	625 x 680 x 580
External dimensions WxHxD (mm)	734 x 905 x 810
Internal Volume (l)	240
Temp. variation across chamber ($\pm^{\circ}\text{C}$)	2.5
Temp. variation at a point ($\pm^{\circ}\text{C}$)	0.3
Number of Shelves	2
Number of Shelf Positions	8
Maximum load per shelf (kg)	20

Model: 279D

Voltage V/hz	230/50
Power W	3000
Shipping Size WxHxD (mm)	920 x 1170 x 1030
Mass (kg)	132
Internal dimensions WxHxD (mm)	725 x 780 x 755
External dimensions WxHxD (mm)	834 x 1005 x 985
Internal Volume (l)	400
Temp. variation across chamber ($\pm^{\circ}\text{C}$)	2.5
Temp. variation at a point ($\pm^{\circ}\text{C}$)	0.3
Number of Shelves	2
Number of Shelf Positions	10
Maximum load per shelf (kg)	20

Model: 276

Voltage V/hz	230/50
Power W	1000
Shipping Size WxHxD (mm)	890 x 600 x 760
Mass (kg)	45
Internal dimensions WxHxD (mm)	485 x 400 x 455
External dimensions WxHxD (mm)	750 x 540 x 690
Internal Volume (l)	80
Temp. variation across chamber ($\pm^{\circ}\text{C}$)	1.5
Temp. variation at a point ($\pm^{\circ}\text{C}$)	0.3
Number of Shelves	2
Number of Shelf Positions	5
Maximum load per shelf (kg)	26

Model: 278

Voltage V/hz	230/50
Power W	2000
Shipping Size WxHxD (mm)	1050 x 980 x 920
Mass (kg)	88
Internal dimensions WxHxD (mm)	625 x 680 x 580
External dimensions WxHxD (mm)	890 x 820 x 810
Internal Volume (l)	240
Temp. variation across chamber ($\pm^{\circ}\text{C}$)	2.5
Temp. variation at a point ($\pm^{\circ}\text{C}$)	0.3
Number of Shelves	2
Number of Shelf Positions	8
Maximum load per shelf (kg)	20

Model: 279

Voltage V/hz	230/50
Power W	2000
Shipping Size WxHxD (mm)	1130 x 1120 x 1030
Mass (kg)	132
Internal dimensions WxHxD (mm)	725 x 780 x 755
External dimensions WxHxD (mm)	990 x 920 x 985
Internal Volume (l)	400
Temp. variation across chamber ($\pm^{\circ}\text{C}$)	2.5
Temp. variation at a point ($\pm^{\circ}\text{C}$)	0.3
Number of Shelves	2
Number of Shelf Positions	10
Maximum load per shelf (kg)	20