



# HEALTHCARE TECHNOLOGIES

## icanCLAVE

Order Code	ICANCLAVE C208	ICANCLAVE C212	ICANCLAVE C218	ICANCLAVE C223
Model Number	STE-8-C	STE-12-C	STE-18-C	STE-23-C
Chamber Dimension	Ø170x320 mm	Ø200x360 mm	Ø247x350 mm	Ø247x450 mm
Chamber Volume	8 L	12 L	18 L	23 L
Number of Trays	2/ Max. 3	2/ Max. 3	3/ Max. 5	3/ Max. 5
Material of Trays	Aluminum	Stainless Steel	Aluminum	Aluminum
Size of Trays.	130x265mm	156x295mm	192x285mm	192x380mm
Voltage Freq.	220/110V 50/60Hz	220/110V 50/60Hz	220/110V 50/60Hz	220/110V 50/60Hz
POWER.r	2000W	2000W	2000W	2200W
Overall Dimension (WxHxD) Weight	420x370x525 mm 37 kgs	420x370x595 mm 40 kgs	490x455x600 mm 47 kgs	490x455x690 mm 53 kgs



HEALTHCARE  
TECHNOLOGIES

icanCLAVE

### CYCLE INFORMATION

Program	Temp. [·C]	Pressure [Mpa]	Sterilization time [min.]
SOLID	134	210	4
	121	110	20
LIQUID	134	210	10
	121	110	30
WRAPPED	134	210	10
	121	110	30
TEXTILE	134	210	10
	121	110	30
PRION	134	210	18
B&D TEST	134	210	3.5
HELIX TEST	134	210	3.5

### STANDARDS

- 93/42/EEC Directive on Medical Devices.
- 97/23/EEC Directive on Pressure Equipment EN 13060 Small Steam Sterilizers.
- EN 61326, EN 61000 Electromagnetic Compatibility.
- EN 61010-1, UL 61010-1 Safety requirements.
- EN 61010-2-040, IEC 61010-2-040. Specific requirements for steam sterilizers used for treatment of medical materials.