

## All-In-One Ristocetin CoFactor Assay System

**vW Select<sup>®</sup>** is comprised of optimized reagents, control plasmas, and lot specific technical information for use in the detection and evaluation of von Willebrand Disease (vWD), including Type 2B and the Quebec Platelet Disorder.



- Quantify Ristocetin CoFactor Activity
- Achieve More Linear Standard Curves
- Monitor Patient Response to Therapy
- Improve Test Performance & Repeatability
- Reduce Assay & Laboratory Variables
- Meet CLIA & Other Regulatory Requirements
- Save Time & Money

**vW Select<sup>®</sup>** delivers superior repeatability in preparing standard reference curves as well as patient testing when compared to other commercial or Lab Developed Tests.

**COMPLETE TEST KIT**

#### NET CONTENTS:

Lyophilized Platelets	2 x 10.0mL
AggRecetin <sup>®</sup>	2 x 15.0mg
vW Normal Reference Plasma	5 x 0.5mL
vW Abnormal Control Plasma	5 x 0.5mL
vW Normal Control Plasma	5 x 0.5mL
Tris Buffered Saline	3 x 10.0mL

An ISO 13485:2003 Registered Company

**Worldwide:** +1 215-441-4000  
**USA:** 1-800-257-3282  
**Fax Worldwide:** +1 215-443-8820  
**customer.service@biodatacorp.com**