

Impactis M

shockwave therapy unit



SHOCKWAVE THERAPY UNIT

- dedicated transmitter for aesthetic medicine
- applicator with an integrated spring shock absorber absorbing vibrations
- ergonomic silhouette of the transmitter cap improving holding comfort
- the device in its original version allows for over 2 000 000 shocks
- the applicator regeneration kit allows for another 2 000 000 or more shocks

SHOCKWAVE EMISSION MODES

- single
- continuous
- burst

TREATMENT PROGRAMS

- 43 preset treatment programs



TECHNICAL PARAMETERS

- | | |
|---|----------------------------------|
| • compressor pressure | 1 – 5 bar |
| • impulse emission frequency (shocks) | 1 – 25 Hz |
| • shocks | 1 – 10 000 |
| • power density | max 3 mJ/mm ² |
| • projectile system life | > 2 million cycles |
| • supply | 100-240 V, 47-63 Hz, 24VDC 6,25A |
| • device dimensions | 36,1 x30,4 x15,1 cm |
| • controller weight | 7 kg |
| • weight with the applicator and power unit | 9 kg |



STANDARD ACCESSORIES

- mains cable (1 pc.)
- switch mode power supply (1 pc.)
- shockwave applicator (1 pc.)
- transmitter 10 mm (1 pc.)
- transmitter 15 mm (1 pc.)
- transmitter 20 mm (1 pc.)
- spare gaskets 8x3, 12x3, 13x3 (2 pcs. each)
- spare elastomer spring (2 pcs.)
- ultrasound gel 500 g (1 pc.)
- spare time-lag fuse (1 pc.)
- applicator holder (1 pc.)
- projectile chamber spanner (1 pc.)
- screwdriver for fixing the holder (1 pc.)
- display cloth (1 pc.)
- display stylus (1 pc.)
- user guide (1 pc.)

OPTIONAL ACCESSORIES

- transmitter 35 mm
- regeneration kit (projectile system)
- hearing protector
- bag for the unit and accessories
- ELF X trolley